

## REFERENCE GUIDE


### PRODUCT PICTURE




OLED PANEL

POWER BUTTON

### DISPLAY




### REAR PANEL



## QUICK START GUIDE

**\*NOTE:**  
WITH THE POWER ON, QUICK PRESSES OF THE POWER BUTTON WILL TOGGLE THE SCREEN DISPLAY. A LONGER PRESS OF THE BUTTON WILL ACCESS THE SETTINGS MENU. FOR FULL INSTRUCTIONS, PLEASE REFER TO THE ENCLOSED USER GUIDE.


01



### LOAD THE BATTERIES

OPEN THE BATTERY COVER BY SLIDING IN THE DIRECTION OF THE GREEN ARROW. INSERT THE 2 AAA BATTERIES, ENSURING THE + AND - ARE THE CORRECT WAY UP AS PER THE MARKINGS ON THE DEVICE. THEN CLOSE THE BATTERY COVER BY SLIDING IT BACK AGAIN IN THE OPPOSITE DIRECTION.

02



### OPEN THE JAWS

SQUEEZE ON THE TOP AND BOTTOM OF THE DEVICE AT THE END WITH THE AVRO LOGO AND PIVOT PIN, TO OPEN THE OXIMETER'S JAWS LIKE A CROCODILE'S MOUTH.

03



### INSERT FINGER INTO OXIMETER

PUSH THE MIDDLE FINGER ON YOUR OTHER HAND (IT DOESN'T MATTER WHICH HAND, BUT MOST USE THE LEFT FOR MEASURING) INTO THE JAWS OF THE DEVICE. DON'T FORCE IT, BUT IT NEEDS TO BE IN FAR ENOUGH SO AS TO BE IN LINE WITH THE PIVOT PIN. YOUR FINGER CAN BE EITHER WAY UP. WITH THE OXIMETER ON YOUR FINGER, REST YOUR ARM ON A FLAT SURFACE AT AROUND CHEST HEIGHT.

04



### OBSERVE YOUR READINGS

PRESS THE POWER BUTTON. THERE WILL BE A 'BEEP' AND YOU WILL SEE A QUICK 'OXIMETER' MESSAGE AND THE DATE. IF IT'S BEEN SET. WITHIN 5-10 SECONDS THE DEVICE WILL SETTLE AND GIVE YOU ACCURATE SPO2 AND BPM READINGS\*. ONCE YOU REMOVE YOUR FINGERTIP, THE DEVICE WILL SHUTDOWN AUTOMATICALLY WITHIN A FEW SECONDS.

\* UNUSUALLY LOW SPO2 READINGS OR HIGH/LOW BPM READINGS WILL BE ACCOMPANIED BY A REGULAR WARNING BEEP.

### \*PLEASE NOTE:

WHILST THIS DEVICE IS EXTREMELY ACCURATE, A COLD OR DAMAGED FINGER WILL SOMETIMES NOT GENERATE AN ACCURATE READING. IF IN DOUBT, TRY YOUR OTHER MIDDLE FINGER, FOLLOWED BY FINGERS OTHER THAN YOUR MIDDLE ONES UNTIL YOU ARE SATISFIED THAT THE READINGS ARE CONSISTENT.